


*EmbraceAirKing*

*EmbraceAirPlus*

*LumbAirPlus*

*InlineSeating*

**Instruction Manual  
& Warranty Guide**



**YOUR NEW INNOTECH SUPPORT PRODUCT**

Congratulations on your wise purchase of Innotech's high quality hand crafted support product. You have invested in a product that is manufactured in Canada from the highest quality raw materials.


By following the enclosed instructions your Innotech support will last many, many years.

**HOW TO USE YOUR INNOTECH BACK SUPPORT**

1. Place the support behind you in a position that feels comfortable.
2. Secure it to the seat with the use of the adjustable buckle and strap system.
3. To inflate, squeeze the bulb repeatedly until the desired firmness is obtained. To release air simply depress the push button valve. Do not exceed 30 pumps.



**USING OUR PATENTED VERTICALLY ADJUSTABLE AIR SYSTEM**

**NOTE:** Always release air completely before you pull either adjustment strap.




**ATTACHING OPTIONAL INLINE SEATING**

**Step 1**  
Attach the enclosed velcro strips inside the cover on rear side directly on the wooden frame of the EmbraceAirPlus, EmbraceAirKing or EmbracePlus.



**Step 2**  
Locate velcro straps on InlineSeating.



**Step 3**  
Evenly attach InlineSeating straps to velcro on EmbraceAirPlus, EmbraceAirKing or EmbracePlus. Press firmly.

**FREQUENTLY ASKED QUESTIONS AND TROUBLESHOOTING**

**Why do I need the Patented Adjustment Feature?**  
Each body is different. Furthermore, when you sit in a car, on a kitchen chair or couch, the seat pan will either sit you up (hard) or allow you to sink (soft). This changes the height of the position of the natural curve in your lumbar spine. Our patented vertical adjustment helps you to compensate for these changes and precisely locate your lumbar spine with these variations in height. Remember always to adjust the air chamber when it is empty.

**How many times do I squeeze the bulb?**  
Until you feel comfortable, the more you squeeze the bulb, the firmer the support will become. Thirty pumps is quite firm. Do not exceed 30 pumps.

**Can I adjust the tube length?**  
Yes, the back support comes to you with ample length for ease of operation. The tube is long enough to sit comfortably on your lap. If you wish to shorten the tube simply remove the bulb and valve, cut at the desired length and reattach the bulb and valve.

**Is there latex in the product?**  
Yes, a small amount is used in the fabrication of the air bladder.

**How long will the air stay in?**  
The firmness is determined by you. Changes can occur through temperature fluctuations, excessive pressure changes or simply through time. Slight changes are natural but generally the firmness should remain constant for hours. It is recommended, dependent upon personal preferences, that the firmness be altered occasionally. This usually changes the micro circulation patterns and may alleviate discomfort or muscle stiffness due to prolonged or extended sitting.

**CLEANING AND MAINTENANCE**


Our exclusive fabric, *Innotex™*, is available only on Innotech's back supports. It has been treated to become immune to bacteria, fungus and odour. It wicks moisture away better than silk, cotton, wool or synthetics. To dust, simply vacuum. For spot cleaning, use any recommended high quality upholstery cleaner. Be sure to follow the manufacturers directions.



Innotech Rehabilitation Products Inc.  
P.O. Box 534,  
Orillia, Ontario, Canada  
L3V 6K2

**Innotech Guarantee**

Your hand crafted individually inspected back support is guaranteed to be free of defects in material and workmanship for a period of 90 days from date of purchase. A **Life Time Guarantee** on the air valve components is offered. In order to be eligible for the guarantee the attached warranty card must be completed and mailed within 15 days after date of purchase. Warranty repairs are **free** but all shipping charges are the responsibility of the customer. For enquiries please call 1-800-361-0228 Eastern Standard Time.



**Important**

To activate your warranty, fill in the following information and mail this card back to us. Warranty must be postmarked within 15 days of purchase.

**Warranty Validation**

Last Name and Initial .....  
Street Address .....  
City ..... Prov./State .....  
Postal Code/Zip ..... Phone (.....).....  
**Purchase Information**  
Date of Purchase Year ..... Month ..... Day .....  
Where did you purchase the product?  
Name .....  
Address .....  
City ..... Prov./State .....  
**Where did you hear about Innotech's support products?**  
Physical Therapist ..... Chiropractor .....  
M.D. .... Pharmacy ..... Other Store .....  
Ads ..... Word of mouth ..... Other .....  
Is this your first Innotech support product? Yes ..... No .....  
If "No", how many others have you owned? .....  
Would you like information about our other products? .....

**Innotech Warranty**

Innotech guarantees this support to be free from defects in material and workmanship for a period of 90 days from date of purchase.

If you find a defect, please contact Innotech directly at 1-800-361-0228.

**Repair Procedures**

- 1 - Contact Innotech, Orillia, Ontario, Canada.
- 2 - Upon receipt of authorization, return the product as directed, including proof and date of purchase, at your expense and risk. Repairs under warranty will be made at our expense.
- 3 - Product will be returned to you at your cost and risk.

**Limitations**

This warranty is void if the product is damaged by improper or abnormal use or by accident.